



PHE-921 pH/ATC Waterproof Combination Electrode

Operator's Manual: M2893/0398



WARRANTY

OMEGA only warrants that the parts manufactured by it will be as specified and free of defects upon receipt by purchaser. OMEGA MAKES NO OTHER WARRANTY OR REPRESENTATIONS OF ANY KIND WHATSOEVER, EXPRESSED OR IMPLIED, EXCEPT THAT OF TITLE AND ALL IMPLIED WARRANTIES INCLUDING ANY WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE HEREBY DISCLAIMED. THE REMEDIES OF PURCHASER SET FORTH HEREIN SHALL NOT EXCEED THE PURCHASE PRICE OF THE PRODUCT. IN NO EVENT SHALL OMEGA BE LIABLE FOR CONSEQUENTIAL, INCIDENTAL, OR SPECIAL DAMAGES.

RETURN REQUESTS / INQUIRIES

Direct all warranty and repair requests/inquiries to the OMEGA ENGINEERING Customer Service Department. BEFORE RETURNING ANY PRODUCT(S) TO OMEGA, PURCHASER MUST OBTAIN AN AUTHORIZED RETURN (AR) NUMBER FROM OMEGA'S CUSTOMER SERVICE DEPARTMENT (IN ORDER TO AVOID PROCESSING DELAYS). The assigned AR number should then be marked on the outside of the return package and on any correspondence.

FOR WARRANTY RETURNS, please have the following information available BEFORE contacting OMEGA:

1. P.O. number under which the product was PURCHASED,
2. Model and serial number of the product under warranty, and
3. Repair instructions and/or specific problems relative to the product.

FOR NON-WARRANTY REPAIRS OR CALIBRATION, consult OMEGA for current repair/calibration charges. Have the following information available BEFORE contacting OMEGA:

1. P.O. number to cover the COST of the repair/calibration,
2. Model and serial number of product, and
3. Repair instructions and/or specific problems relative to the product.

OMEGA's policy is to make running changes, not model changes, whenever an improvement is possible. This affords our customers the latest in technology and engineering. OMEGA is a registered trademark of OMEGA ENGINEERING, INC. © Copyright 1995 OMEGA ENGINEERING, INC. All rights reserved. This documentation may not be copied, photocopied, reproduced, translated, or reduced to any electronic medium or machine-readable form, in whole or in part, without prior written consent of OMEGA ENGINEERING, INC.

Germany/Austria:
Dammierstrasse 26, D-75392
Deckenfronm, Germany
Tel: 49 (07056) 3017 FAX: 49 (07056) 8540
Toll Free in Germany: 0130 11 21 66
e-mail: germany@omega.com

United Kingdom: ISO 9002 Certified
25 Swannington Road, Broughton Ashtey,
Leicestershire, LE9 6TL, England
Tel: 44 (1455) 285520 FAX: 44 (1455) 283912
P.O. Box 7, Omega Drive, Htlam,
Manchester, M44 5EX, England
Tel: 44 (161) 777-6611 FAX: 44 (161) 777-6622
Toll Free in England: 0800-488-488
e-mail: uk@omega.com

Benelux:
Postbus 8034, 1180 LA Amstelveen,
The Netherlands
Tel: (31) 20 6418405 FAX: (31) 20 6434643
Toll Free in Benelux: 06 0993344
e-mail: nl@omega.com

Czech Republic:
ul. Rude armady 1868, 733 01 Karvina-Hranice,
Czech Republic
Tel: 420 (69) 6311899 FAX: 420 (69) 6311114
e-mail: czech@omega.com

France:
9, rue Denis Papin, 78190 Trappes
Tel: (33) 130-621-400 FAX: (33) 130-699-120
Toll Free in France: 0800-4-06342
e-mail: france@omega.com

Servicing Europe:

USA: ISO 9001 Certified
One Omega Drive, Box 4047
Stamford, CT 06907-0047
Tel: (203) 359-1660
FAX: (203) 359-7700
e-mail: info@omega.com

Servicing North America:


Canada:
976 Bergar
Laval (Quebec) H7L 5A1
Tel: (514) 856-6928
FAX: (514) 856-6886
e-mail: canada@omega.com

Mexico and Latin America:
Tel: (95) 800-TC-OMEGA
FAX: (95) 203-359-7807
Customer Service: 1-800-622-BEST
Engineering Service: 1-800-872-9436 /
1-800-USA-WHEN
EASLINK: 62968934 CABLE: OMEGA

For immediate technical or application assistance:

Mexico and Latin America:
Tel: (95) 800-TC-OMEGA
FAX: (95) 203-359-7807
E-mail: espanol@omega.com

OMEGA^{net} On-Line Service
<http://www.omega.com>
Internet e-mail
info@omega.com





PHE-921 pH/ATC Waterproof Combination Electrode

Operator's Manual: M2893/0398

Specifications

Range	0 - 14 pH
Isopotential Point	0 ± 30 mV
Temperature Range	0 to 100°C
Reference Element	Silver/Silver Chloride
Junction	Ceramic
Filling Solution	4M KCl saturated with AgCl
Storage Solution	4M KCl
Electrode Diameter	12mm
Body Length	108mm
Cable Length	1 meter
Connector Type	BNC and 2.5mm phone plug
Response Time	95% in 10 seconds
Sodium Error	<0.2 pH at pH 14 with 1M sodium

Electrode Preparation

1. Remove the electrode from Storage Vial.
2. Clean any salt deposits from the electrode body by rinsing with deionized (DI) or distilled water.
3. Verify that the Electrode Filling Solution is at least 2.5" (64mm) high. Must keep Filling Solution above the Fill Mark. If not, slide the electrode Fill Cap down to uncover the Fill Holes, and using the plastic dropper fill the inner chamber with Filling Solution.
4. Connect the electrode to the meter as described in the meter's manual.
5. During use leave the Fill Cap in the open position.

Measuring Techniques

- Standardize at least daily for the most accurate readings.
- Always use fresh buffers.
- Stir all buffers and samples.
- During standardization, allow time for the electrode to stabilize before entering the buffer into the meter.
- Rinse the electrode with DI water between samples.
- Blot the electrode dry (do NOT rub or wipe).

Storage

Place the electrode in its Storage Vial containing Electrode Storage Solution. Electrode Storage Solution (4M KCl) can be purchased or made. Make 4M KCl solution by dissolving 7.45 grams of potassium chloride into 100 mL of DI or distilled water. Do NOT store the electrode in deionized water.

Troubleshooting

1. Connect the pH electrode to the BNC connector.
2. Place the meter into absolute (*not* relative) mV mode.
3. Rinse and place the clean electrode into a pH 7 buffer:
The mV reading should be 0 ± 30 mV.
Outside of this range indicates a problem.
See below for solutions.*
4. Rinse and place the clean electrode into a pH 4 buffer:
The mV reading should be 159 to 186 mV more than value in pH 7.
Example: pH 7 buffer is 10 mV; pH 4 buffer should read 169 to 196 mV.
Outside of this range indicates a problem.
See below for solutions.*
5. Optionally, rinse and place the clean electrode into a pH 10 buffer:
The mV reading should be -159 to -186 mV less than value in pH 7.
Example: pH 7 buffer is 10 mV; pH 10 buffer should read -149 to -176 mV.
Outside of this range indicates a problem.
See below for solutions.*

* Ensure that buffers are fresh and not contaminated. If readings are still outside of this range, clean the electrode and repeat test.

Cleaning Procedures

If the electrode becomes sluggish, or has been determined to be bad, clean the electrode by following the directions below.

General

Soak in 0.1M HCl or 0.1M HNO₃ for half an hour.
Drain and refill the Electrode Filling Solution.
Soak the electrode in Filling or Storage Solution for 1 hour.

Inorganic

Soak in 0.1M tetrasodium EDTA solution 15 minutes.
Drain and refill the Electrode Filling Solution.
Soak the electrode in Filling or Storage Solution for 1 hour.

Protein

Soak in 1% pepsin in 0.1M HCl for 15 minutes.
Drain and refill the Electrode Filling Solution.
Soak the electrode in Filling or Storage Solution for 1 hour.

Grease and Oil

Rinse with mild detergent or methanol solution.
Drain and refill the Electrode Filling Solution.
Soak the electrode in Filling or Storage Solution for 1 hour.