

# **ISE-8755 and ISE 8756**

## **Silver/Sulfide Ion Electrodes**



**Operator's Manual**



## Servicing USA and Canada: Call OMEGA Toll Free

### USA

One Omega Drive, Box 4047  
Stamford, CT 06907-0047  
Telephone: (203) 359-1660  
FAX: (203) 359-7700

### Canada

976 Bergar  
Laval (Quebec) H7L 5A1  
Telephone: (514) 856-6928  
FAX: (514) 856-6886

Sales Service: 1-800-826-6342 / 1-800-TC-OMEGA<sup>SM</sup>

Customer Service: 1-800-622-2378 / 1-800-622-BEST<sup>SM</sup>

Engineering Service: 1-800-872-9436 / 1-800-USA-WHEN<sup>SM</sup>

TELEX: 996404 EASYLINK: 62968934 CABLE: OMEGA

OMEGAnet<sup>SM</sup> On-Line Service

<http://www.omega.com>

## Servicing Europe: United Kingdom Sales and Distribution Center

25 Swannington Road, Broughton Astley, Leicestershire, LE9 6TU, England

Telephone: 44 (1455) 285520 or 0800-488-488 FAX: 44 (1455) 283912

## The OMEGA Complete Measurement and Control Handbooks & Encyclopedias

- ✓ Temperature
- ✓ Pressure, Strain & Force
- ✓ Flow and Level
- ✓ pH and Conductivity
- ✓ Data Acquisition Systems
- ✓ Electric Heaters
- ✓ Environmental Monitoring and Control



**Call for Your FREE Handbook Request Form Today: (203) 359-RUSH**



Remove the Packing List and verify that you have received all equipment, including the following (quantities in parentheses):

- ISE-8755 or ISE-8756 Electrode (1)
- Sheet of semilog paper (1)
- Sheet of fine polishing paper (1)
- Operator's Manual (1)

If you have any questions about the shipment, please call the OMEGA Customer Service Department. When you receive the shipment, inspect the container and equipment for signs of damage. Note any evidence of rough handling in transit. Immediately report any damage to the shipping agent.

---

**NOTE**

The carrier will not honor damage claims unless all shipping material is saved for inspection. After examining and removing contents, save packing material and carton in the event reshipment is necessary.

---

	Page
<b>Section 1 General Instructions .....</b>	<b>1-1</b>
1.1 Introduction .....	1-1
1.2 Required Equipment .....	1-1
1.3 Required Solutions .....	1-1
1.3.1 For Silver .....	1-1
1.3.2 For Sulfide .....	1-2
<b>Section 2 General Preparation .....</b>	<b>2-1</b>
2.1 Electrode Preparation .....	2-1
2.2 Electrode Slope Check With a Standard pH/mV Meter .....	2-1
2.2.1 Using Silver Standards .....	2-1
2.2.2 Using Sulfide Standards .....	2-1
2.3 Electrode Slope Check With an Ion Meter .....	2-1
2.3.1 Using Silver Standards .....	2-2
2.3.2 Using Sulfide Standards .....	2-2
<b>Section 3 Measurement .....</b>	<b>3-1</b>
3.1 Measuring Hints .....	3-1
3.2 Sample Requirements .....	3-1
3.3 Units of Measurement .....	3-2
<b>Section 4 Measurement Procedure .....</b>	<b>4-1</b>
4.1 Direct Measurement .....	4-1
4.1.1 Direct Measurement of Silver Using a Standard pH/mV Meter .....	4-1
4.1.2 Direct Measurement of Silver Using an Ion Meter .....	4-2
4.1.3 Direct Measurement of Sulfide Using a Standard pH/mV Meter .....	4-3
4.1.4 Direct Measurement of Sulfide Using an Ion Meter .....	4-4
4.2 Low-Level Silver Determinations Using a pH/mV Meter .....	4-5
4.3 Titration .....	4-6
4.3.1 Titration of Sulfide .....	4-6
4.3.2 Titration of Silver .....	4-7
4.3.3 Low-Level Chloride Titration .....	4-8
4.3.4 Indicator Titration .....	4-10

	<b>Page</b>
<b>Section 5 Electrode Characteristics .....</b>	<b>5-1</b>
5.1 Reproducibility .....	5-1
5.2 Interferences .....	5-1
5.3 Complexation .....	5-1
5.4 Temperature Influences .....	5-1
5.5 Electrode Response .....	5-2
5.6 Limits of Detection .....	5-3
5.7 pH Effects .....	5-4
5.8 Electrode Life .....	5-4
5.9 Electrode Storage .....	5-4
<b>Section 6 Electrode Theory and Operation .....</b>	<b>6-1</b>
<b>Section 7 Troubleshooting Guide .....</b>	<b>7-1</b>
7.1 Meter .....	7-1
7.2 Glassware .....	7-1
7.3 Electrodes .....	7-1
7.4 Standards and Reagents .....	7-2
7.5 Sample .....	7-2
7.6 Technique .....	7-2
7.7 Troubleshooting Hints .....	7-2
<b>Section 8 Specifications and Ordering Information .....</b>	<b>8-1</b>

**Notes**

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

## 1.1 Introduction

OMEGA® Silver/Sulfide Ion Electrodes are used to quickly, simply, accurately, and economically measure silver and sulfide ions in aqueous solutions. The two ions are virtually never present in solution together, owing to the extreme insolubility of silver sulfide. Indirect measurements of cyanide or measurements of halide ions by titration may be done with this electrode.

## 1.2 Required Equipment

1. A pH/mV meter or an ion meter, either line-operated or portable.
2. Semi-logarithmic, four-cycle graph paper for preparing calibration curves when using the meter in the mV mode (included with the electrode). Gran's plot paper (10% volume corrected) is recommended for low-level chloride measurement.
3. A magnetic stirrer.
4. The OMEGA® Silver/Sulfide Ion Electrode, ISE-8755 (reference electrode necessary), or the OMEGA Silver/Sulfide Ion Combination Glass Electrode, ISE-8756. Use OMEGA Filling Solution (Part Number PHFS-B250) in the ISE-8756.
5. The OMEGA® Double Junction Reference Electrode, PHE-3211 (for use with the ISE-8755), with OMEGA Filling Solution, PHFS-A250, in the inner junction and OMEGA Filling Solution, PHFS-B250, in the outer junction.
6. Polishing Paper to polish dirty or etched electrode membranes (one piece supplied with the electrode).

## 1.3 Required Solutions

Deionized or distilled water for solution and standard preparation. Water used in the preparation of sulfide standards and of SAOB should also be deaerated.

### 1.3.1 For Silver

1. OMEGA's Ionic Strength Adjuster (ISA), Part Number ISE-8755-R1. To prepare this solution from your own laboratory stock, half fill a 1000-mL volumetric flask with distilled water and add 425 grams of reagent-grade sodium nitrate,  $\text{NaNO}_3$ . Swirl the flask to dissolve the solid. Fill the flask to the mark with distilled water, cap, and upend the flask several times to mix the contents.
2. OMEGA's Silver Standard Solution, 0.1M  $\text{AgNO}_3$ , Part Number ISE-8755-S1. To prepare this solution from your own laboratory stock, dry reagent-grade, pulverized silver nitrate in a laboratory oven for one hour at 150°C. Quantitatively transfer 16.99 grams of the dried silver nitrate to a one-liter flask containing about 500 mL of distilled water. Swirl the flask to dissolve the solid. Fill the flask to the mark with distilled water, cap, and upend several times to mix the solution. Store the solution in a brown bottle, tightly capped, in a dark place.
3. OMEGA's Silver Standard, 1000 ppm as  $\text{Ag}^+$ , Part Number ISE-8755-S2. To prepare this solution from your own laboratory stock, dry reagent-grade, pulverized



silver nitrate in a laboratory oven for one hour at 150°C. Quantitatively transfer 1.57 grams of the dried silver nitrate to a one-liter volumetric flask containing about 500 mL of distilled water. Swirl the flask to dissolve the solid. Fill the flask to the mark with distilled water, cap, and upend several times to mix the solution. Store the solution in a brown bottle, tightly capped, in a dark place.

4. Chloride Titrant (low-level),  $2.82 \times 10^{-3}$  M  $\text{AgNO}_3$ . Dry reagent-grade, pulverized silver nitrate as described above. Quantitatively transfer 0.479 gram of dried silver nitrate to a one-liter volumetric flask about half-full of distilled water. Swirl the flask to dissolve the solid. Fill to the mark with distilled water, cap, and upend several times to mix the solution. Store in a brown bottle in a dark place.

### 1.3.2 For Sulfide

1. Sulfide Anti-Oxidant Buffer (SAOB). This buffer must be used fresh and will range in color from clear to yellow-brown. It has become oxidized when it turns dark brown and should then be discarded. Store fresh SAOB in a tightly stoppered bottle. To prepare SAOB from your own stock, fill a one-liter volumetric flask with about 500 mL of distilled, deaerated water, 200 mL of 10 M NaOH, 35 grams of ascorbic acid, and 67 grams of disodium EDTA. Swirl the mixture until the solids dissolve and fill to the mark with distilled, deaerated water.
2. Lead Perchlorate Solution, 0.1M. This solution is used for the titration of sulfide standards.
3. Sulfide Standard. To prepare a stock solution of saturated sodium sulfide, add about 100 grams of reagent grade  $\text{Na}_2\text{S} \cdot 9\text{H}_2\text{O}$  to 100 mL of distilled, deaerated water. Shake well and let stand overnight, storing in a tightly stoppered bottle in a hood.

Prepare a weekly sulfide standard by adding 500 mL SAOB to a one-liter volumetric flask, pipetting 10 mL of the stock solution into the flask and diluting to the mark with distilled, deaerated water.

The exact concentration (C) can be determined by titrating 10 mL of the standard with 0.1 M lead perchlorate. Use the silver/sulfide ion electrode (and the reference electrode) to indicate the endpoint.

The calculation is as follows:

$$C = 3206 (V_t/V_s),$$

where:

C = concentration as ppm sulfide

$V_t$  = volume of titrant at endpoint

$V_s$  = volume of standard used (10 mL).

Prepare other standards **each day** by serial dilution of the weekly standard. To do a tenfold dilution, accurately measure 10 mL of the standard and add it to a 100-mL volumetric flask. Add 45 mL of SAOB and dilute to the mark with distilled, deaerated water.



## 2.1 Electrode Preparation

Remove the rubber caps covering the electrode tips and the rubber insert covering the filling hole of the reference electrode. Fill the combination electrode or the reference electrode with the filling solution shipped with the electrode to a level just below the fill hole. No preparation is required for a sealed reference electrode. Connect the electrodes to the proper terminals as recommended by the meter manufacturer.

## 2.2 Electrode Slope Check With a pH/mV Meter

Check electrodes each day. Slope is defined as the change in potential observed when the concentration changes by a factor of 10.

### 2.2.1 Using Silver Standards

1. To a 150-mL beaker, add 100 mL of distilled water and 2 mL of ISA. Place the beaker on a magnetic stirrer and begin stirring at a constant rate. After ensuring that the meter is in the mV mode, lower the electrode tips into the solution.
2. Using a pipet, add 1 mL of 0.1 M or 1000-ppm silver standard to the beaker. When the reading has stabilized, record the millivolt reading.
3. Using a pipet, add 10 mL of the same silver standard used above to the beaker. When the reading has stabilized, record the millivolt reading.
4. Determine the difference between the two readings. A difference of  $57 \pm 2$  mV indicates correct electrode operation, assuming the solution temperature is between 20 and 25°C. Refer to the Troubleshooting section if the potential change is not within this range.

### 2.2.2 Using Sulfide Stock Solutions

1. To a 150-mL beaker, add 50 mL of distilled water and 50 mL of SAOB solution. Place the beaker on a magnetic stirrer and stir at a constant rate. After ensuring that the meter is in mV mode, lower the electrode tips into the solution.
2. Using a pipet, add 1 mL of weekly standard sulfide solution to the beaker. When the reading has stabilized, record the millivolt reading.
3. Using a pipet, add 10 mL of weekly standard sulfide used above to the beaker. When the reading has stabilized, record the millivolt reading.
4. Determine the difference between the two readings. A difference of  $26 \pm 2$  mV indicates correct electrode operation, assuming the solution temperature is between 20 and 25°C. Refer to the Troubleshooting section if the potential change is not within this range.

## 2.3 Electrode Slope Check With an Ion Meter

Check electrodes each day.

### 2.3.1 Using Silver Standards:

1. Prepare standard silver solutions with concentrations that vary by tenfold. Use either the 0.1 M  $\text{AgNO}_3$  or the 1000-ppm stock solutions. Use the serial dilution method for this preparation.
2. To a 150-mL beaker, add 100 mL of the lower value standard and 2 mL of ISA. Place the beaker on a magnetic stirrer and begin stirring at a constant rate. Lower the electrode tips into the solution. Ensure that the meter is in the concentration mode.
3. Adjust the meter to the concentration of the standard and fix the value in the memory according to the meter manufacturer's instructions.
4. Rinse the electrodes with distilled water and blot dry.
5. To another 150-mL beaker, add 100 mL of the higher value standard and 2 mL of ISA. Place the beaker on a magnetic stirrer and begin stirring at a constant rate. Lower the electrode tips into the solution.
6. Adjust the meter to the concentration of the standard and fix the value in memory.
7. Read the electrode slope according to the meter manufacturer's instructions. Correct electrode operation is indicated by a slope of 90–100%. Refer to the Troubleshooting section if the potential is not within this range.

### 2.3.2 Using Sulfide Stock Solutions

1. Prepare standard sulfide solutions with concentrations that vary by tenfold. Use the serial dilution method and the weekly sulfide standard to prepare these solutions.
2. To a 150-mL beaker, add 50 mL of the lower value standard, 25 mL of SAOB, and 25 mL of distilled water. Place the beaker on a magnetic stirrer and begin stirring at a constant rate. Lower the electrode tips into the solution. Ensure that the meter is in the concentration mode.
3. Adjust the meter to the concentration of the standard and fix the value in the memory according to the meter manufacturer's instructions.
4. Rinse the electrodes with distilled water and blot dry.
5. To a 150-mL beaker, add 50 mL of the higher value standard, 25 mL of SAOB, and 25 mL of distilled water. Place the beaker on a magnetic stirrer and begin stirring at a constant rate. Lower the electrode tips into the solution.
6. Adjust the meter to the concentration of the standard and fix the value in memory.
7. Read the electrode slope according to the meter manufacturer's instructions. Correct electrode operation is indicated by a slope of 90–100%. Refer to the Troubleshooting section if the slope is not within this range.

### 3.1 Measuring Hints

All samples and standards should be at the same temperature for precise measurement. A difference of one degree C in temperature will result in about a 2% measurement error.

Constant, but not violent, stirring is necessary for accurate measurement. Magnetic stirrers can generate sufficient heat to change the solution temperature. To counteract this effect, place a piece of insulating material, such as a styrofoam sheet, between the stirrer and beaker.

All silver samples and silver standards should be stored away from light. Always dilute sulfide samples with SAOB (1:1) when they are collected.

Always rinse the electrodes with distilled water and blot dry between measurements. Use a clean, dry tissue to prevent cross-contamination.

A slow-responding electrode may point to deposits on the membrane. Use polishing paper to remove such deposits. Refer to the section on Electrode Response.

For low level silver measurements, use labware of plastic. For samples with high ionic strength, prepare standards with a composition similar to the sample. Always check to see that the membrane is free from air bubbles after immersion into the standard or sample.

### 3.2 Sample Requirements

All samples must be aqueous and must not contain organics that can dissolve the epoxy electrode body and/or the cement that bonds the sensing crystal to the electrode body. Infrequent measurements in solutions containing methanol, benzene, or acetone are permitted. Please check with OMEGA before using the electrode in other organic solvents.

The temperature of the standard solution and of the sample solution should be the same and should be below 80°C. Samples containing sulfide must be buffered with SAOB to convert  $\text{HS}^{-1}$  and  $\text{H}_2\text{S}$  to  $\text{S}^{-2}$ . Samples containing silver must be below pH 8 to avoid reaction with  $\text{OH}^{-1}$ . Silver samples should be acidified with 1 M  $\text{HNO}_3$ , if necessary.

Mercury cannot be present in silver samples. Since  $\text{HgS}$  and  $\text{Hg}_2\text{S}$  are insoluble, no mercury will be present in sulfide samples. Other interferences should be absent. If they are present, use the procedures found in the Interferences and Electrode Response sections to remove them.

### 3.3 Units of Measurement

Silver and sulfide concentrations are measured in units of parts per million, equivalents per liter, moles per liter, or any other convenient concentration unit. Table 1 indicates some of the concentration units.

**Table 1. Concentration Unit Conversion Factors**

ppm $S^{-2}$	ppm $Ag^{+1}$	N( $S^{-2}$ )	M( $Ag^{+1}$ )
321.0	1079.0	$2.00 \times 10^{-2}$	$1.00 \times 10^{-2}$
100.0	—	$6.23 \times 10^{-3}$	$3.12 \times 10^{-3}$
32.1	108.0	$2.00 \times 10^{-3}$	$1.00 \times 10^{-3}$
—	100.0	$1.84 \times 10^{-3}$	$9.27 \times 10^{-4}$
3.2	10.8	$2.00 \times 10^{-4}$	$1.00 \times 10^{-4}$

## 4.1 Direct Measurement

Direct measurement is a simple procedure for measuring a large number of samples. A single meter reading is all that is required for each sample. The ionic strength of samples and standards should be made the same by adjustment with SAOB for all sulfide solutions and with ISA for all silver solutions. The temperature of both the sample solution and standard solution should be the same.

### 4.1.1 Direct Measurement of Silver Using a pH/mV Meter

1. By serial dilution of the 0.1 M or 1000-ppm standards, prepare  $10^{-2}$ ,  $10^{-3}$ , and  $10^{-4}$  M or 100- and 10-ppm silver standards. Add 2 mL of ISA per 100 mL of standard. Prepare standards with a composition similar to that of the samples, if the samples have an ionic strength above 0.1 M.
2. Place the most dilute solution ( $1.0 \times 10^{-4}$  M or 10 ppm) on the magnetic stirrer and begin stirring at a constant rate. After ensuring that the meter is in the mV mode, lower the electrode tips into the solution. When the reading has stabilized, record the mV reading.
3. Place the midrange solution ( $1.0 \times 10^{-3}$  M or 100 ppm) on the magnetic stirrer and begin stirring. After rinsing the electrodes with distilled water, blot dry and immerse the electrode tips in the solution. When the reading has stabilized, record the mV reading.
4. Place the most concentrated solution ( $1.0 \times 10^{-2}$  M or 1000 ppm) on the magnetic stirrer and begin stirring. After rinsing the electrodes with distilled water, blot dry and immerse the electrode tips in the solution. When the reading has stabilized, record the mV reading.
5. Using the semi-logarithmic graph paper, plot the mV reading (linear axis) against concentration (log axis). Extrapolate the curve down to about  $2.0 \times 10^{-6}$  M or 0.2 ppm. A typical calibration curve can be found in Figure 1.

A calibration curve is constructed on semi-logarithmic paper when using a pH/mV meter in the millivolt mode. The measured electrode potential in mV (linear axis) is plotted against the standard concentration (log axis). In the linear region of the curve, only three standards are necessary to determine a calibration curve. In the non-linear region, additional points must be measured. The direct measurement procedures given are for the linear portion of the curve. The non-linear portion of the curve requires the use of low-level procedures.

6. To a clean, dry 150-mL beaker, add 100 mL of sample and 2 mL of ISA. Place the beaker on the magnetic stirrer and begin stirring. Place the electrode tips in the solution. When the reading has stabilized, record the mV reading. Determine the concentration directly from the calibration curve.




Figure 1. Typical Silver Electrode Calibration Curve

7. The calibration should be checked every two hours. Assuming no change in ambient temperature, place the electrode tips in the midrange standard. After the reading has stabilized, compare it to the original reading recorded in Step 3 above. A reading differing by more than 0.5 mV or a change of ambient temperature will necessitate the repetition of Steps 2–5 above. A new calibration curve should be prepared daily.

#### 4.1.2 Direct Measurement of Silver Using an Ion Meter

1. By serial dilution of the 0.1 M or 1000-ppm silver standard, prepare two silver standards with a concentration near the expected sample concentration. Measure out 100 mL of each standard into individual 150-mL beakers and add 2 mL of ISA to each.
2. Place the more dilute solution on the magnetic stirrer and begin stirring at a constant rate. Ensure that the meter is in the concentration mode.
3. Lower the electrode tips into the solution.
4. Adjust the meter to the concentration of the silver standard and fix the value in the memory according to the meter manufacturer's instructions after stabilization of the reading.
5. Rinse the electrodes with distilled water and blot dry.
6. Place the more concentrated solution on the magnetic stirrer and begin stirring at a constant rate.
7. Lower the electrode tips into the solution.

8. Adjust the meter to the concentration of the silver standard and fix the value in memory according to the meter manufacturer's instructions after stabilization of the reading.
9. For low-level measurements, place the rinsed, dried electrodes into the solution containing 100 mL of distilled water and 2 mL of ISA. After stabilization, fix the blank value in the meter according to the meter manufacturer's instructions.
10. Place 100 mL of the sample and 2 mL of ISA in a 150-mL beaker. Place it on the magnetic stirrer and begin stirring.
11. Immerse the electrode tips in the solution and wait for the reading to stabilize. Read the concentration directly from the meter display.
12. The calibration should be checked every two hours. Assuming no change in ambient temperature, place the electrode tips in the first silver standard. After the reading has stabilized, compare it to the original reading in Step 4 above. A reading differing by more than 0.5 mV or a change in the ambient temperature will necessitate the repetition of Steps 2–8(9) above. The meter should be recalibrated daily.

#### 4.1.3 Direct Measurement of Sulfide Using a pH/mV Meter

1. By serial dilution of the weekly standard, prepare three standard sulfide solutions. Measure out 50 mL of each standard into individual 150-mL beakers and add 25 mL of SAOB and 25 mL of distilled water to each.
2. Place the most dilute solution on the magnetic stirrer and begin stirring at a constant rate. After ensuring that the meter is in the mV mode, lower the electrode tips into the solution. When the reading has stabilized, record the mV reading.
3. Place the midrange sulfide standard on the magnetic stirrer and begin stirring. After rinsing the electrodes with distilled water, blot dry and immerse the electrode tips in the solution. When the reading has stabilized, record the mV reading.
4. Place the most concentrated sulfide standard on the magnetic stirrer and begin stirring. After rinsing the electrodes with distilled water, blot dry and immerse the electrode tips in the solution. When the reading has stabilized, record the mV reading.
5. Using the semi-logarithmic graph paper, plot the mV reading (linear axis) against concentration (log axis). A typical calibration curve can be found in Figure 2.

A calibration curve is constructed on semilogarithmic paper when using a pH/mV meter in the millivolt mode. The measured electrode potential in mV (linear axis) is plotted against the standard concentration (log axis). In the linear region of the curve, only three standards are necessary to determine a calibration curve. In the non-linear region, additional points must be measured. The direct measurement procedures given are for the linear portion of the curve. The non-linear portion of the curve requires the use of low-level procedures.




Figure 2. Typical Sulfide Electrode Calibration Curve

6. To a clean, dry 150-mL beaker, add 50 mL of the sulfide sample, 25 mL of SAOB, and 25 mL of distilled water. Place the beaker on the magnetic stirrer and begin stirring. Place the electrode tips in the solution. When the reading has stabilized, record the mV reading. Determine the concentration directly from the calibration curve.
7. The calibration should be checked every two hours. Assuming no change in ambient temperature, place the electrode tips in the midrange standard. After the reading has stabilized, compare it to the original reading recorded in Step 3 above. A reading differing by more than 0.5 mV or a change in the ambient temperature will necessitate the repetition of Steps 2–5 above. A new calibration curve should be prepared daily.

#### 4.1.4 Direct Measurement of Sulfide Using an Ion Meter

1. By serial dilution of the weekly standard, prepare two sulfide standards with a concentration near the expected sample concentration. Measure out 50 mL of each standard into individual 150-mL beakers and add 25 mL of SAOB and 25 mL of distilled water to each.
2. Place the more dilute solution on the magnetic stirrer and begin stirring at a constant rate. Ensure that the meter is in the concentration mode.
3. Lower the electrode tips into the solution.
4. Adjust the meter to the concentration of the sulfide standard and fix the value in memory according to the meter manufacturer's instructions after stabilization of the reading.



5. Rinse the electrodes with distilled water and blot dry.
6. Place the more concentrated solution on the magnetic stirrer and begin stirring at a constant rate.
7. Lower the electrode tips into the solution.
8. Adjust the meter to the concentration of the sulfide standard and fix the value in the memory according to the meter manufacturer's instructions after stabilization of the reading.
9. For low-level measurements, place the rinsed, dried electrodes into the solution containing 50 mL of SAOB and 50 mL of distilled water. After stabilization, fix the blank value in the meter according to the meter manufacturer's instructions.
10. Place 50 mL of the sample, 25 mL of SAOB, and 25 mL of distilled water in a 150-mL beaker, place it on the magnetic stirrer, and begin stirring.
11. Immerse the electrode tips in the solution and wait for the reading to stabilize. Read the concentration directly from the meter display.
12. The calibration should be checked every two hours. Assuming no change in ambient temperature, place the electrode tips in the first sulfide standard. After the reading has stabilized, compare it to the original reading in Step 4 above. A reading differing by more than 0.5 mV or a change in the ambient temperature will necessitate the repetition of Steps 2–8(9) above. The meter should be recalibrated daily.

## 4.2 Low-Level Silver Determinations Using a pH/mV Meter

This procedure is recommended for solutions with ionic strengths less than  $1.0 \times 10^{-2}$  M. If the solution is high in ionic strength, but low in silver, use the same procedure, but prepare a calibration solution with a composition similar to the samples.

1. Using 1 mL of standard ISA, dilute to 100 mL with distilled water. This low-level ISA ( $1.0 \text{ M NaNO}_3$ ) is added at the rate of 1 mL low-level ISA to each 100 mL of sample. The background ionic strength will be  $1.0 \times 10^{-2}$  M.
2. Dilute 1 mL of 0.1 M standard to one liter to prepare a  $1.0 \times 10^{-4}$  M solution for measurements in moles per liter. Prepare a 10-ppm standard solution by diluting 1 mL of the 1000-ppm standard to 100 mL of solution for measurements in ppm. Standards should be prepared fresh daily.
3. Add 1 mL of low-level ISA to a 100-mL volumetric flask and fill to the mark with distilled water. Pour this solution into a 150-mL beaker and place the beaker on the magnetic stirrer. Begin stirring at a constant rate.
4. Place the electrode tips in the solution. Ensure that the meter is in the mV mode.
5. Add increments of the  $1.0 \times 10^{-4}$  M or 10-ppm silver standard as given in Table 2 below.
6. After the reading has stabilized, record the mV reading.

**Table 2. Stepwise Calibration for Low-Level Silver Measurements**

Step	Pipet	Added Volume (mL)	Concentration ppm	M
1	A	0.1	0.01	$1.0 \times 10^{-7}$
2	A	0.1	0.02	$2.0 \times 10^{-7}$
3	A	0.2	0.04	$4.0 \times 10^{-7}$
4	A	0.2	0.06	$6.0 \times 10^{-7}$
5	A	0.4	0.10	$9.9 \times 10^{-7}$
6	B	2.0	0.29	$2.9 \times 10^{-6}$
7	B	2.0	0.48	$4.8 \times 10^{-6}$

Pipet A = 1-mL graduated pipet

Pipet B = 2-mL pipet

Solutions: Additions of 10-ppm or  $1.0 \times 10^{-4}$  M standard to 100 mL of ISA, as prepared in Step 3 above

7. On semi-logarithmic graph paper, plot the millivolt reading (linear axis) against the concentration (log axis) as in Figure 2.
8. Rinse the electrodes in distilled water and blot dry.
9. Measure out 100 mL of the sample into a 150-mL beaker, add 1 mL of low-level ISA, and place the beaker on the magnetic stirrer. Begin stirring. Lower the electrode tips into the solution. After the reading has stabilized, record the mV reading and determine the concentration from the low-level calibration curve.
10. Prepare a new low-level calibration curve daily. Check the calibration curve every two hours by repeating Steps 2–7.

## 4.3 Titration

Titration is a very accurate determination of total sulfide or total silver concentration. This method makes use of the electrode as an endpoint detector. The electrode can also be used to determine halide concentrations.

### 4.3.1 Titration of Sulfide

The minimum sulfide sample concentration for this method is  $1.0 \times 10^{-5}$  M. The titrant to be used is a lead perchlorate standard solution.

1. Using Lead Perchlorate Solution, 0.1 M, prepare a lead titrant that is about 10–20 times as concentrated as the expected sample concentration by dilution.
2. Dilute 25 mL of the sample with 25 mL of SAOB in a 150-mL beaker. Place the beaker on the magnetic stirrer and begin stirring. Lower the tips of the electrodes into the solution.

3. Using a 10-mL burette, add titrant in 0.5- to 1.0-mL increments. Record the mV reading against the volume of titrant added. As the mV potential change increases, add smaller increments, down to 0.1- to 0.2-mL increments. Continue to add titrant and record the mV potential against the volume until little change is noted in the mV reading, even when adding 0.5- to 1.0-mL increments.
4. Using linear graph paper, plot the mV readings (y-axis) against the volume (x-axis). The endpoint is determined at the steepest slope on the titration curve, as illustrated in Figure 3.




Figure 3. Typical Titration of 25 mL of 0.03 M Sulfide with 0.1 M  $\text{Pb}(\text{ClO}_4)_2$

5. The sample concentration,  $C_s$ , is calculated before the dilution with SAOB, as follows:

$$C_s = (V_t/V_s) C_t,$$

where:

$V_s$  = sample volume before dilution (25 mL)

$V_t$  = titrant volume at endpoint

$C_t$  = titrant concentration (M).

#### 4.3.2 Titration of Silver

The Silver/Sulfide Ion Electrode is a highly sensitive endpoint detector for silver titrations with a halide standard solution. It can also be used as an indicator for halide titrations with a silver standard solution.

The minimum silver sample concentration for this method is about  $1.0 \times 10^{-4}$  M. For halides, it is  $1.0 \times 10^{-4}$  M.

Data gathered by titration of less concentrated samples can be plotted using the Gran's Plot technique. (See the chloride titration that follows.)

1. Using Sodium Chloride Solution, 0.1 M NaCl, for silver titrations, or the 0.1 M silver standard for halide titrations, prepare a titrant about 10–20 times the expected sample concentration by dilution.
2. Measure out 50 mL of the sample into a 150-mL beaker, place the beaker on the magnetic stirrer, and begin stirring. Place the tips of the electrodes in the solution.
3. Using a 10-mL burette, add titrant in 0.5- to 1.0-mL increments. Record the mV reading against the volume of titrant added. As the mV potential change increases, add smaller increments, down to 0.1- to 0.2-mL increments. Continue to add titrant and record the mV potential against the volume until little change is noted in the mV reading, even when adding 0.5- to 1.0-mL increments.
4. Using linear graph paper, plot the mV readings (y-axis) against the volume (x-axis). The endpoint is determined at the steepest slope on the titration curve as illustrated in Figure 3.
5. The sample concentration ( $C_s$ ) is calculated as follows:

$$C_s = (V_t/V_s) C_t,$$

where:

$V_s$  = sample volume (50 mL)

$V_t$  = titrant volume at endpoint

$C_t$  = titrant concentration (M).

### 4.3.3 Low-Level Chloride Titration

Chloride samples with a concentration below  $1.0 \times 10^{-4}$  M or 5 ppm cannot be titrated in the usual manner since the solubility product of silver chloride is not exceeded.

Using 10% Gran's Plot paper, data from samples down to  $4.0 \times 10^{-6}$  (0.04 ppm) may be plotted.

1. Using 20 mL of standard ISA, dilute to 100 mL with distilled water. This low-level ISA (1.0 M NaNO<sub>3</sub>) is added at the rate of 1 mL low-level ISA to each 100 mL of sample. The background ionic strength will be  $1.0 \times 10^{-2}$  M.
2. Measure out 100 mL of distilled water into a 150-mL beaker, add 1 mL of low-level ISA, place the beaker on the magnetic stirrer, and begin stirring.
3. Prepare a 10-mL burette with the low-level chloride titrant ( $2.82 \times 10^{-3}$  M AgNO<sub>3</sub>) mentioned in the General Instructions (Required Solutions) section.
4. Lower the tips of the electrodes into the solution.
5. Add 1 mL of titrant and record the mV reading after stabilization. Repeat this operation four more times.

6. Using Gran's Plot paper, plot the values as follows:
  - a. The horizontal axis is the volume axis. Each major division is equal to 1 mL of titrant used.
  - b. The vertical axis is the mV axis. Each major division is equal to 5 mV. The axis should be scaled so that the mV value obtained after 5 mL of titrant added is near the top of the graph.
7. After drawing a straight line through the points, the line should intersect the horizontal axis at 0 mL. (See Figure 4.)
8. Measure out 100 mL of the sample into a 150-mL beaker, add 1 mL of low-level ISA, place the beaker on the magnetic stirrer, and begin stirring.
9. Prepare a 10-mL burette with the low-level chloride titrant mentioned in Step 3 above.
10. Allow the titrant to run slowly into the sample until the mV reading is near the value recorded for the first 1-mL increment added in Step 5 above. Record the mV reading and the total volume added.
11. Add 1-mL increments of titrant and record the mV reading vs the total volume until the mV reading approaches the highest mV reading recorded in Step 5 above.
12. Plot the mV reading vs total volume for all points on the Gran's Plot paper. Draw a straight line through as many points as possible, extending the line to the horizontal axis. The volume intercept value is equal to the chloride concentration in ppm.




Figure 4. Gran's Plot Titration of Chloride Sample With  $\text{AgNO}_3$ , Using a Silver Electrode

#### 4.3.4 Indicator Titration

A special method for measuring low-level cyanide concentration. The silver / sulfide electrode is capable of cyanide measurements down to 0.03 ppm. (Midgley, D.; Torrance, K., 1991, *Potentiometric Water Analysis*, 2nd Edition, p. 406–413, John Wiley & Sons, Chichester.)

## 5.1 Reproducibility

Electrode measurements reproducible to 2% can be obtained if the electrode is calibrated every hour. Factors such as temperature fluctuations, drift, and noise limit reproducibility. Reproducibility is independent of concentration within the electrode's operating range.

## 5.2 Interferences

A surface layer of silver metal may be formed by strongly reducing solutions. A layer of silver salt may be deposited on the membrane if high levels of ions forming very insoluble salts are present in the sample. Performance may be restored by polishing. See the section on Electrode Response for the proper polishing procedure.

All silver samples must be free of mercury. Sulfide samples will not have mercury present due to the extreme insolubility of HgS and Hg<sub>2</sub>S. Biological samples and protein in food interfere with silver measurements, but the protein interference can be removed by acidifying to pH 2–3 with 1 M HNO<sub>3</sub>.

## 5.3 Complexation

Total concentration ( $C_t$ ) whether sulfide or silver ions, consists of free ions ( $C_f$ ) and complexed or bound ions ( $C_c$ ) in solution:

$$C_t = C_f + C_c.$$

Since the electrode only responds to free ions, any complexing agent in the solution reduces the measured concentration of ions.

Silver ions complex with many species, notably cyanide, thiosulfate, ammonia, and chelants, such as EDTA.

Sulfide ions form complexes with hydrogen ions (HS<sup>-1</sup> and H<sub>2</sub>S). Sulfide ions also form soluble complexes with elemental sulfur (S<sub>2</sub><sup>-2</sup>, S<sub>3</sub><sup>-2</sup>, S<sub>4</sub><sup>-2</sup>, etc.) and tin, arsenic, and antimony ions.

## 5.4 Temperature Influences

Samples and standards should be within ±1°C of each other, since electrode potentials are influenced by changes in temperature. Because of the solubility equilibria on which the electrode depends, the absolute potential of the reference electrode changes slowly with temperature. The slope of the electrode, as indicated by the factor "S" in the Nernst equation, also varies with temperature. Table 3 gives values for the "S" factor in the Nernst equation for each ion.

**Table 3. Temperature vs. Values for the Electrode Slope**

Temp (°C)	S <sup>-2</sup>	Ag <sup>+1</sup>
0	27 • 1	54 • 2
10	28 • 1	56 • 2
20	29 • 1	58 • 2
25	29 • 6	59 • 2
30	30 • 1	60 • 1
40	31 • 1	62 • 1
50	32 • 1	64 • 1

If changes in temperature occur, the electrodes should be recalibrated.

The temperature range for OMEGA's Silver/Sulfide Ion Electrode is 0–80°C, provided that temperature equilibrium has occurred. If the temperature varies substantially from room temperature, equilibrium times up to one hour are recommended.

## 5.5 Electrode Response

Plotting the electrode mV potential against the silver concentration on semi-logarithmic paper results in a straight line with a slope of about 57 mV per decade. (Refer to Figure 1, page 4-2.) The sulfide ion also gives a straight line when the electrode mV potential is plotted against the sulfide concentration, but the slope is about 26 mV per decade. (See Figure 2, page 4-4.)

The time needed to reach 99% of the stable electrode potential reading (the electrode response time) varies from several seconds in highly concentrated solutions to several minutes near the detection limit. (Refer to Figure 5.)

A drifting potential reading or a decrease in electrode slope may mean that the electrode membrane needs polishing.

To polish the membrane:

1. If using polishing paper, cut off a one- to two-inch piece and place it face up on the laboratory bench.
2. Put a few drops of distilled or deionized water in the center of the paper.
3. Holding the paper (cotton) steady with one hand, bring the membrane of the electrode down perpendicular to the paper, and with a slight swirling motion, gently polish the tip of the electrode against the surface of the polishing paper (cotton) for a few seconds.
4. Rinse the electrode surface with distilled or deionized water and soak the electrode tip in standard solution for about five minutes before use.






Figure 5. Typical Electrode Time Response to Step Changes in  $\text{AgNO}_3$

5. If using jeweller's rouge, place a cotton ball on the table top and flatten it using the bottom of a beaker.
6. Put one to two drops of distilled or deionized water in the center of the cotton pad.
7. Add a small amount of jeweller's rouge to the damp cotton.
8. Continue with Steps 3 and 4 above.

## 5.6 Limits of Detection

The upper limit of detection in pure silver nitrate solutions is 1 M. In the presence of other ions, the upper limit of detection is above  $1.0 \times 10^{-1}$  M silver, but two factors influence this upper limit. Both the possibility of a liquid junction potential developing at the reference electrode and the salt extraction effect influence this upper limit. Some salts may extract into the electrode membrane at high salt concentrations, causing deviation from the theoretical response. Either dilute samples between 1 M and  $1.0 \times 10^{-1}$  M or calibrate the electrode at four or five intermediate points.

The lower limit of detection is influenced by the slight water solubility of the electrode pellet. Refer to Figures 1 and 2 for a comparison of the theoretical response to the actual response at low levels of sulfide and silver.

## 5.7 pH Effects

Silver reacts with hydroxide ions to form a precipitate of  $\text{Ag}_2\text{O}$ , in ammonia-free basic solutions. By keeping all solutions slightly acidic, this can be avoided. Adjust the pH of silver solutions below 8, if necessary, with 1 M  $\text{HNO}_3$ .

Bisulfide ion ( $\text{HS}^-$ ) and hydrogen sulfide ( $\text{H}_2\text{S}$ ) result when hydrogen ion complexes sulfide ion. Larger amounts of sulfide ion are complexed as the pH is lowered. The use of SAOB in all samples containing sulfide maintains a fixed level of  $\text{S}^{2-}$  ions, since the free sulfide ion ( $\text{S}^{2-}$ ) exists in only very basic solutions. In the acid range, sulfide is chiefly in the form of  $\text{H}_2\text{S}$ , while in the pH range 6–12, almost all the sulfide is in the  $\text{HS}^-$  form.

## 5.8 Electrode Life

The silver/sulfide electrode will last six months in normal laboratory use. On-line measurements might shorten operational lifetime to several months. In time, the response time will increase and the calibration slope will decrease to the point calibration is difficult and electrode replacement is required.

## 5.9 Electrode Storage

OMEGA's Silver/Sulfide Electrodes may be stored for short periods of time in  $1.0 \times 10^{-2}$  M silver (or sulfide) solution. For longer storage (more than two weeks), rinse and dry the sensing pellet and cover the membrane tip with any protective cap shipped with the electrode. The reference portion of the combination electrode (or the outer chamber of the reference electrode) should be drained of filling solution, if refillable, and the rubber insert should be placed over the filling hole.

The OMEGA Silver/Sulfide Ion Electrode is composed of a silver sulfide crystal membrane bonded into an epoxy or glass body. When an electrode potential develops across the membrane, the membrane is in contact with a solution containing sulfide or silver ions and is capable of measuring free sulfide or silver ions. This electrode potential is measured against a constant reference potential, using a pH/mV meter or an ion meter. The level of sulfide or silver ions, corresponding to the measured potential, is described by the Nernst equation:

$$E = E_o + S \log X,$$

where:

$E$  = measured electrode potential

$E_o$  = reference potential (a constant)

$S$  = electrode slope (-26 mV for sulfide; +57 mV for silver)

$X$  = level of sulfide or silver in solution.

The activity ( $X$ ) represents the effective concentration of the ions in solution. The activity is related to the free ion concentration ( $C_f$ ) by the activity coefficient ( $\gamma$ ) by:

$$X = \gamma C_f.$$

Activity coefficients vary, depending on total ionic strength ( $I$ ), defined as:

$$I = 1/2 \sum C_x Z_x^2,$$

where:

$C_x$  = concentration of ion  $X$

$Z_x$  = charge of ion  $X$

$\Sigma$  = sum of all of the types of ions in the solution.


In the case of high and constant ionic strength relative to the sensed ion concentration, the activity coefficient ( $\gamma$ ) is constant and the activity ( $X$ ) is directly proportional to the concentration.

To adjust the background ionic strength to a high and constant value, ionic strength adjuster (ISA) is added to samples and standards. The recommended solution for sulfide is SAOB, used to prevent oxidation and free sulfide ion from hydrogen ion, in addition to adjusting the ionic strength. The recommended ISA for silver is  $\text{NaNO}_3$ . Solutions other than these may be used as ionic strength adjusters as long as ions that they contain do not interfere with the electrode's response to sulfide ions or to silver ions. Samples with high ionic strength (greater than 0.1 M) do not need ISA added and standards for these solutions should be prepared with a composition similar to the samples.

The reference electrode must also be considered. When two solutions of different composition are brought into contact with one another, liquid junction potentials arise. Millivolt potentials occur from the interdiffusion of ions in the two solutions. Electrode charge will be carried unequally across the solution boundary resulting in a potential difference between the two solutions, since ions diffuse at different rates. When making measurements, it is important to remember that this potential be the same when the reference is in the standardizing solution, as well as in the sample solution or the change in liquid junction potential will appear as an error in the measured electrode potential.

The composition of the liquid junction filling solution in the reference electrode is most important. The speed with which the positive and negative ions in the filling solutions diffuse into the sample should be as nearly equal as possible, that is, the filling solution should be equitransferent. No junction potential can result if the rate at which positive and negative charge carried into the sample solution is equal.

Strongly acidic ( $\text{pH} = 0\text{--}2$ ) and strongly basic ( $\text{pH} = 12\text{--}14$ ) solutions are particularly troublesome to measure. The high mobility of hydrogen and hydroxide ions in samples make it impossible to mask their effect on the junction potential with any concentration of an equitransferent salt. One must either calibrate the electrodes in the same pH range as the sample or use a known increment method for ion measurement.



# Troubleshooting Guide

---

The goal of troubleshooting is the isolation of a problem through the checking of each of the system components in turn: the meter, the glassware, the electrodes, the standards and reagents, the sample, and the technique.

## 7.1 Meter

The meter may be checked by following the checkout procedure in the instrument instruction manual.

## 7.2 Glassware

Clean glassware is essential for good measurement. Be sure to wash the glassware well with a mild detergent and rinse very well with distilled or deionized water. Clean glassware will drain without leaving water droplets behind.

## 7.3 Electrodes

The electrodes may be checked by using the procedure found in the Electrode Slope Check section.

1. Be sure to use distilled or deionized water when following the procedures given in Electrode Slope Check.
2. If the electrode fails to respond as expected, see the sections Measuring Hints and Electrode Response. Repeat the slope check.
3. If the electrodes still fail to respond as expected, substitute another silver/sulfide ion electrode that is known to be in good working order for the questionable electrode. If the problem persists and you are using an electrode pair, try the same routine with a working reference electrode.
4. If the problem persists, the reagent may be of poor quality, interferences in the sample may be present, or the technique may be faulty. (Refer to the Standards and Reagents, Sample, and Technique sections below.)
5. If another electrode is not available for test purposes, or if the electrode in use is suspect, review the instruction manual and be sure to:
  - Clean and rinse the electrodes thoroughly
  - Prepare the electrodes properly
  - Use the proper filling solution
  - Adjust the pH and the ionic strength of the solution by the use of the proper ISA
  - Measure correctly and accurately.
  - Review the Troubleshooting Hints section.

## 7.4 Standards and Reagents

Whenever problems arise with the measuring procedure that has been used successfully in the past, be sure to check the standard and reagent solutions. If in doubt about the credibility of any of the solutions, prepare them again. Errors may result from contamination of the ISA, incorrect dilution of standards, poor quality distilled/deionized water, or a simple mathematical miscalculation.

## 7.5 Sample

Look for possible interferences, complexing agents, or substances that could affect the response or physically damage the sensing electrode (or the reference electrode) if the electrodes work perfectly in the standard, but not in the sample.

Try to determine the composition of the samples prior to testing to eliminate a problem before it starts. (Refer to the Measuring Hints, Sample Requirements, and Interferences sections.)

## 7.6 Technique

Be sure that the electrode's limit of detection has not been exceeded. Be sure that the analysis method is clearly understood and is compatible with the sample.

Refer to the instruction manual again. Reread the General Preparation and Electrode Characteristics sections.

If trouble still persists, call the pH Department of OMEGA Engineering.

## 7.7 Troubleshooting Hints

Symptom	Possible Causes	Next Step
<b>Out of Range Reading</b>	Defective meter	Check meter with shorting strap (see meter instruction manual)
	Defective electrode	Check electrode operation
	Electrode not plugged in properly	Unplug electrode and reseal
	Reference electrode not filled	Be sure that reference electrode is filled
	Air bubble on membrane	Remove bubble by redipping electrode
	Electrodes not in solution	Put electrodes in solution

Symptom	Possible Causes	Next Step
<b>Noisy or Unstable Readings (readings continuously or randomly changing)</b>	Defective meter	Check meter with shorting strap
	Air bubble on membrane	Remove bubble by redipping electrode
	ISA not used	Use recommended ISA
	Meter or stirrer not grounded	Ground meter or stirrer
	Defective electrode	Replace electrode
	Electrode exposed to interferences	Soak electrode in silver or sulfide standard
<b>Drift (reading slowly changing in one direction)</b>	Samples and standards at different temperatures	Allow solutions to come to room temperature before measurement
	Electrode exposed to complexing agents	Check section entitled Complexation
	Incorrect reference filling solution	Use recommended filling solution
<b>Low Slope or No Slope</b>	Standards contaminated or incorrectly made	Prepare fresh standards
	ISA not used	Use recommended ISA
	Standard used as ISA	Use ISA
	Electrode exposed to complexing agents	Check section entitled Complexation
	Air bubble on membrane	Remove bubble by redipping probe

Symptom	Possible Causes	Next Step
<b>"Incorrect Answer" (but calibration curve is good)</b>	Incorrect scaling of semi-log paper	Plot millivolts on the linear axis; on the log axis, be sure concentration numbers within each decade are increasing with increasing concentration
	Incorrect sign	Be sure to note the sign of the millivolt number correctly
	Incorrect standards	Prepare fresh standards
	Wrong units used	Apply correct conversion factor: $1.0 \times 10^{-3} \text{ M} = 32.1 \text{ ppm S}^{-2} =$ $2.0 \times 10^{-3} \text{ N (S}^{-2}\text{)}; 1.0 \times 10^{-3} \text{ M} =$ $108 \text{ ppm Ag}^{+1}$
	Complexing agents in sample	Check entitled Complexation section



## Specifications

<b>Concentration Range:</b>	1 M to $1.0 \times 10^{-7}$ M sulfide (32,000 to 0.003 ppm) 1 M to $1.0 \times 10^{-7}$ M silver (108,000 to 0.01 ppm)
<b>pH Range:</b>	2 to 12
<b>Temperature Range:</b>	0–80°C
<b>Resistance:</b>	<1 megaohm
<b>Reproducibility:</b>	±2%
<b>Samples:</b>	Aqueous solutions only; no organic solvents
<b>Size:</b>	110-mm length 12-mm diameter 1.0-m cable length
<b>Storage:</b>	Store in silver or sulfide solution

## Ordering Information

Part Number	Description
<b>ISE-8755</b>	Silver/Sulfide Ion Electrode, Mono (reference electrode necessary), Epoxy Body
<b>ISE-8756</b>	Silver/Sulfide Ion Electrode, Combination, Glass Body
<b>PHE-3211</b>	Reference Electrode, Double Junction, Sleeve, Epoxy Body, for Use With the ISE-8755
<b>ISE-8755-S1</b>	Silver Standard, 0.1M AgNO <sub>3</sub>
<b>ISE-8755-S2</b>	Silver Standard, 1000 ppm AgNO <sub>3</sub>
<b>ISE-8755-R1</b>	Silver Ionic Strength Adjuster (ISA), 5 M NaNO <sub>3</sub>
<b>PHFS-A250</b>	Reference Electrode Filling Solution, 4 M KCl/Ag <sup>+</sup> , for the 5731429 Inner Chamber
<b>PHFS-B250</b>	Reference Electrode Filling Solution, 1 M KNO <sub>3</sub> , for the 5731429 Outer Chamber and ISE-8756 Electrode
<b>948201</b>	Replacement 0.3-micron Aluminum Oxide on Plastic Backing Fine Polishing Paper for the Silver/Sulfide Electrodes



## **WARRANTY/DISCLAIMER**

**OMEGA is glad to offer suggestions on the use of our various products. Nevertheless, OMEGA only warrants that the parts manufactured by it will be as specified and free of defects upon receipt by purchaser.**

**OMEGA MAKES NO OTHER WARRANTIES OR REPRESENTATIONS OF ANY KIND WHATSOEVER, EXPRESSED OR IMPLIED, EXCEPT THAT OF TITLE AND ALL IMPLIED WARRANTIES INCLUDING ANY WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE HEREBY DISCLAIMED.**

**LIMITATION OF LIABILITY: The remedies of purchaser set forth herein are exclusive and the total liability of OMEGA with respect to this order, whether based on contract, warranty, negligence, indemnification, strict liability or otherwise, shall not exceed the purchase price of the component upon which liability is based. In no event shall OMEGA be liable for consequential, incidental or special damages.**

Every precaution for accuracy has been taken in the preparation of this manual; however, OMEGA ENGINEERING, INC. neither assumes responsibility for any omissions or errors that may appear nor assumes liability for any damages that result from the use of the products in accordance with the information contained in the manual.

**SPECIAL CONDITION: Should this equipment be used in or with any nuclear installation or activity, purchaser will indemnify OMEGA and hold OMEGA harmless from any liability or damage whatsoever arising out of the use of the equipment in such a manner.**

## **RETURN REQUESTS/INQUIRIES**

**BEFORE RETURNING ANY PRODUCT(S) TO OMEGA, PURCHASER MUST OBTAIN AN AUTHORIZED RETURN (AR) NUMBER FROM OMEGA'S CUSTOMER SERVICE DEPARTMENT (IN ORDER TO AVOID PROCESSING DELAYS). The assigned AR number should then be marked on the outside of the return package and on any correspondence.**

---

OMEGA's policy is to make running changes, not model changes, whenever an improvement is possible. This affords our customers the latest in technology and engineering.

OMEGA is a registered trademark of OMEGA ENGINEERING, INC.

© Copyright 1996 OMEGA ENGINEERING, INC. All rights reserved. This documentation may not be copied, photocopied, reproduced, translated, or reduced to any electronic medium or machine-readable form, in whole or in part, without prior written consent of OMEGA ENGINEERING, INC.

# Where Do I Find Everything I Need for Process Measurement and Control? **OMEGA...Of Course!**

## **TEMPERATURE**

- ☒ Thermocouple, RTD & Thermistor Probes, Connectors, Panels & Assemblies
- ☒ Wire: Thermocouple, RTD & Thermistor
- ☒ Calibrators & Ice Point References
- ☒ Recorders, Controllers & Process Monitors
- ☒ Infrared Pyrometers

## **PRESSURE/STRAIN FORCE**

- ☒ Transducers & Strain Gages
- ☒ Load Cells & Pressure Gauges
- ☒ Displacement Transducers
- ☒ Instrumentation & Accessories

## **FLOW/LEVEL**

- ☒ Rotameters, Gas Mass Flowmeters & Flow Computers
- ☒ Air Velocity Indicators
- ☒ Turbine/Paddlewheel Systems
- ☒ Totalizers & Batch Controllers

## **pH/CONDUCTIVITY**

- ☒ pH Electrodes, Testers & Accessories
- ☒ Benchtop/Laboratory Meters
- ☒ Controllers, Calibrators, Simulators & Pumps
- ☒ Industrial pH & Conductivity Equipment

## **DATA ACQUISITION**

- ☒ Data Acquisition and Engineering Software
- ☒ Communications-Based Acquisition Systems
- ☒ Plug-in Cards for Apple, IBM & Compatibles
- ☒ Datalogging Systems
- ☒ Recorders, Printers & Plotters

## **HEATERS**

- ☒ Heating Cable
- ☒ Cartridge & Strip Heaters
- ☒ Immersion & Band Heaters
- ☒ Flexible Heaters
- ☒ Laboratory Heaters

## **ENVIRONMENTAL MONITORING AND CONTROL**

- ☒ Metering & Control Instrumentation
- ☒ Refractometers
- ☒ Pumps & Tubing
- ☒ Air, Soil & Water Monitors
- ☒ Industrial Water & Wastewater Treatment
- ☒ pH, Conductivity & Dissolved Oxygen Instruments