____ OMEGA[™]

HHC250

RPM Meter INSTRUCTION

SHEET

M5513/0418 Shop online at omega.com e-mail: info@omega.com For latest product manuals: omegamanual.info

APO HOLD MAX-MIN AVG 188888 CE OMEGA HHC250 RPM METER

INTRODUCTION

This instrument is a non-contact rotational speed meter. The meter measures the rotational speed in RPM (Revolutions Per Minute) and RPS (Revolutions Per Second) with reflective markers stuck on the rotating object to be measured.

Using the meter:


Before start measuring, wait till the rotating object come to completely stops, stick the reflective marker(s) on the object to be measured. (The number of markers can be from 1 to 9. Please refer to MAX/MIN Mode for setting the number of markers).

Attention:

Do not expose the sensor to the light source directly or put reflective tape nearby reflective media, like glass or shiny material to ensure the correctness of reading.

When measuring, hold the meter steady and vertical (within 30°) to the surface of reflective markers.

(A () HOLE



The ideal distance to the measurement object is between 10 and 40 cm.

Suggested reflective tape size:

diagram will be formed as:


Measuring distance between 10 and 20cm tape size 1cm*2cm(H*L).

Measuring distance between 20 and 30cm tape size 2cm*3cm(H*L).

Measuring distance between 30 and 40cm tape size 2cm*5cm(H*L).

Attention:

CE




Reflective tape

Light scatting

B: Red light - effective area.

When the measuring distance is diagram more than 30cm, please align the "B area" to the reflective tape to receive better signal/reading.



SPECIFICATIONS

Measuring range: 100.0 to 99999 rpm

1.7 to 1666 rps **Resolution:** 0.1 rpm (100.0 to 19999.9 rpm)

1 rpm (20000 to 99999 rpm)

Accuracy: Stated accuracy at 23°C±5°C, <75% R.H.

Temperature Coefficient: 0.1 times the applicable accuracy specification per °C out of 23°C±5°C. **Operating environment:** 0°C to 50°C at <70%

R.H.

Storage temperature: -20°C to 60°C, 0 to 80% R.H. with battery removed from meter.

Battery: 1.5Vx2pcs AAA size.

Battery Life: 30 hours typical with alkaline battery.

Low battery indication: The "
"
" is displayed when the battery voltage drops below the operating level.

Dimensions: 24.7mm(T) x 50.9mm(W) x 132.9mm(H).

Weight: Approx. 110g including batteries.

OPERATIONS

There are 3 operation modes, namely Measure-

ment Mode, Set Mode and MAX/MIN Mode.

1. Turning on and off the meter:

When power is off, a short push on "O" key turns on the meter and enter Measurement Mode. When power is on, press and hold "O" key for 2 seconds turns the power off.

Right after power on, all the indicators on the LCD display lights up for one second.

2. Measurement Mode:

Pressing Hold "A" key to measure in rpm or rps. Release "▲" key, and stop measurement and will hold the measuring values.

Backlight:

In measurement mode, a short push on "O" key turns on the backlight for dark environment, it goes off automatically after 15 seconds if without further operation.

3. Set Mode:

In set mode, user can turn on/off APO function and select measurement unit rpm or rps. (APO = Auto Power Off, when APO is ON, themeter power off automatically if no operation in

10 minutes.)

When power is off, press and hold "O" key for 2 seconds to enter Set Mode.

In Set Mode, press "Mode" key each time move settable position in following sequence (indicator of settable position blinks.)

- 3.1 Unit rpm/rps: short push on "A" key switches between rpm and rps. Push "Mode" key to select the unit and enter APO ON/OFF setting.
- 3.2 APO ON/OFF: short push on " \blacktriangle " key switches between ON and OFF. Push "Mode" key to select the APO ON/OFF and exit Set Mode and Entering Measurement Mode. (Please refer to the section of Measurement Mode.)
- XNote: To turn off power in Set Mode abandon the current setting and previous settings remain unchanged. If " appeared, the setting value can works till power off but did not save to meter.

4. MAX/MIN Mode:

In Measurement Mode, short push "Mode" key to enter MAX/MIN Mode. Display MAX, MIN, MAX-MIN, AVG value. For the last time the value of measurement.

Push "Mode" key each time the display changes in following sequence:

4.1 Maximum value recorded: with indicator


"MAX".

- 4.2 Minimum value recorded: with indicator "MIN".
- 4.3 Maximum-Minimum value recorded: with indicator "MAX-MIN".
- 4.4 Average of values recorded: with indicator "AVG".
- 4.5 Set the number of reflective marker (NO.1~9). Press "****" key to scroll the number of reflective markers from 1~9. Push "Mode" key to select the reflective markers and exit the MAX/MIN Mode and return to Measurement Mode.

To clear the recorded values:

In MAX/MIN Mode and during viewing MAX, MIN, MAX-MIN or AVG, press and hold "Mode" key for 2 seconds, to clear the recorded data and return to measurement mode.

MAINTENANCE **Installing and Replacing Battery**




A. Battery Cover B. Battery

- 1. Power is supplied by 2pcs 1.5V (AAA SIZE).
- 2. The "
 "
 "
 appears in the display when battery replacement is needed.
- 3. Push the Battery Cover and lift it in the direction as shown in the figure.
- 4. Remove the batteries from battery compartment.
- 5. Replace with 2 new AAA batteries with polarity as indicated on the bottom of Battery Compartment.
- 6. Replace the Battery Cover.

Cleaning

Periodically wipe the case with a damp cloth and detergent, do not use abrasives or solvents.



OK

0.1 rps (1.7 to 333.3 rps) 1 rps (334 to1666 rps) Accuracy: ±0.02% of reading + 1 dgt **Units:** rpm, rps

WARRANTY/DISCLAIMER

OMEGA ENGINEERING, INC. warrants this unit to be free of defects in materials and workmanship for a period of **13 months** from date of purchase. OMEGA's WARRANTY adds an additional one (1) month grace period to the normal **one (1) year product warranty** to cover handling and shipping time. This ensures that OMEGA's customers receive maximum coverage on each product.

If the unit malfunctions, it must be returned to the factory for evaluation. OMEGA's Customer Service Department will issue an Authorized Return (AR) number immediately upon phone or written request. Upon examination by OMEGA, if the unit is found to be defective, it will be repaired or replaced at no charge. OMEGA's WARRANTY does not apply to defects resulting from any action of the purchaser, including but not limited to mishandling, improper interfacing, operation outside of design limits, improper repair, or unauthorized modification. This WARRANTY is VOID if the unit shows evidence of having been tampered with or shows evidence of having been damaged as a result of excessive corrosion; or current, heat, moisture or vibration; improper specification; misapplication; misuse or other operating conditions outside of OMEGA's control. Components in which wear is not warranted, include but are not limited to contact points, fuses, and triacs.

OMEGA is pleased to offer suggestions on the use of its various products. However, OMEGA neither assumes responsibility for any omissions or errors nor assumes liability for any damages that result from the use of its products in accordance with information provided by OMEGA, either verbal or written. OMEGA warrants only that the parts manufactured by the company will be as specified and free of defects. OMEGA MAKES NO OTHER WARRANTIES OR REPRESENTATIONS OF ANY KIND WHATSOEVER, EXPRESSED OR IMPLIED, EXCEPT THAT OF TITLE, AND ALL IMPLIED WARRANTIES INCLUDING ANY WARRANTY OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE ARE HEREBY DISCLAIMED. LIMITATION OF LIABILITY: The remedies of purchaser set forth herein are exclusive, and the total liability of OMEGA with respect to this order, whether based on contract, warranty, negligence, indemnification, strict liability or otherwise, shall not exceed the purchase price of the component upon which liability is based. In no event shall OMEGA be liable for consequential, incidental or special damages.

CONDITIONS: Equipment sold by OMEGA is not intended to be used, nor shall it be used: (1) as a "Basic Component" under 10 CFR 21 (NRC), used in or with any nuclear installation or activity; or (2) in medical applications or used on humans. Should any Product(s) be used in or with any nuclear installation or activity, medical application, used on humans, or misused in any way, OMEGA assumes no responsibility as set forth in our basic WARRANTY/DISCLAIMER language, and, additionally, purchaser will indemnify OMEGA and hold OMEGA harmless from any liability or damage whatsoever arising out of the use of the Product(s) in such a manner.

RETURN REQUESTS/INQUIRIES

Direct all warranty and repair requests/inquiries to the OMEGA Customer Service Department. BEFORE RETURNING ANY PRODUCT(S) TO OMEGA, PURCHASER MUST OBTAIN AN AUTHORIZED RETURN (AR) NUMBER FROM OMEGA'S CUSTOMER SERVICE DEPARTMENT (IN ORDER TO AVOID PROCESSING DELAYS). The assigned AR number should then be marked on the outside of the return package and on any correspondence.

The purchaser is responsible for shipping charges, freight, insurance and proper packaging to prevent breakage in transit.

FOR **WARRANTY** RETURNS, please have the following information available BEFORE contacting OMEGA:

- Purchase Order number under which the product was PURCHASED,
- 2. Model and serial number of the product under warranty, and
- 3. Repair instructions and/or specific problems relative to the product.

FOR **NON-WARRANTY** REPAIRS, consult OMEGA for current repair charges. Have the following information available BEFORE contacting OMEGA:

- Purchase Order number to cover the COST of the repair,
- 2. Model and serial number of the product, and
- 3. Repair instructions and/or specific problems relative to the product.

OMEGA's policy is to make running changes, not model changes, whenever an improvement is possible. This affords our customers the latest in technology and engineering.

OMEGA is a trademark of OMEGA ENGINEERING, INC.

© Copyright 2018 OMEGA ENGINEERING, INC. All rights reserved. This document may not be copied, photocopied, reproduced, translated, or reduced to any electronic medium or machine-readable form, in whole or in part, without the prior written consent of OMEGA ENGINEERING, INC.

CE OMEGA"

omega.com info@omega.com Servicing North America:

U.S.A.

Headquarters:

Omega Engineering, Inc. Toll-Free: 1-800-826-6342 (USA & Canada only) Customer Service: 1-800-622-2378 (USA & Canada only) Engineering Service: 1-800-872-9436 (USA & Canada only) Tel: (203) 359-1660 Fax: (203) 359-7700 e-mail: info@omega.com

For Other Locations Visit omega.com/worldwide

Where Do I Find Everything I Need for Process Measurement and Control? OMEGA...Of Course! Shop online at omega.com

TEMPERATURE

 Thermocouple, RTD & Thermistor Probes, Connectors, Panels & Assemblies
 Wire: Thermocouple, RTD & Thermistor
 Calibrators & Ice Point References
 Recorders, Controllers & Process Monitors
 Infrared Pyrometers

PRESSURE, STRAIN AND FORCE

Transducers & Strain Gages
 Load Cells & Pressure Gages
 Displacement Transducers
 Instrumentation & Accessories

FLOW/LEVEL

Rotameters, Gas Mass Flowmeters & Flow Computers
 Air Velocity Indicators
 Turbine/Paddlewheel Systems
 Totalizers & Batch Controllers

pH/CONDUCTIVITY

pH Electrodes, Testers & Accessories
 Benchtop/Laboratory Meters
 Controllers, Calibrators, Simulators & Pumps
 Industrial pH & Conductivity Equipment

DATA ACQUISITION

Communications-Based Acquisition Systems
 Data Logging Systems
 Wireless Sensors, Transmitters, & Receivers
 Signal Conditioners
 Data Acquisition Software

HEATERS

Heating Cable
 Cartridge & Strip Heaters
 Immersion & Band Heaters
 Flexible Heaters
 Laboratory Heaters

ENVIRONMENTAL MONITORING AND CONTROL

Metering & Control Instrumentation
 Refractometers
 Pumps & Tubing
 Air, Soil & Water Monitors
 Industrial Water & Wastewater Treatment
 pH, Conductivity & Dissolved Oxygen Instruments